

Shared simplicity

CS 7200 Neo Edition



© 2023 Carestream Dental LLC

CS 7200 Neo Edition

CONFIDENT DIAGNOSIS

Excellent image quality with nearly no risk of over- or underexposure.

FOR EVERYDAY NEEDS AND EVERY PRACTICE

For daily intraoral indications in single or multiple room practices

SCAN & GO SIMPLICITY Enhanced with Scan & Go technology and Smart Plates

AS EASY TO POSITION AS FILM

Flexible, film-like plates with no cable

POWERED WITH CS IMAGING VERSION 8

Provides you with one-stop access to all your images





Comfortable for you

Easy to learn and use


- Easy to install and integrate
- Minimal learning curve
- Film-like positioning

Eliminates darkroom and chemicals

- Faster workflow
- Cost and timesaving
- Safer for the environment

Smarter, more effective workflow

- Scan plates in any order
- No risk of image mix-up





Comfortable for your patient

Safe and fast

- Minimized risk of cross contamination
- Hygienic bags
- Touchless operation

No waiting time

Short reading time

Improved communication

Images are easier to understand

Reduced radiation dose


Lower radiation as compared to film





CS 7200 Neo Edition

AT A GLANCE








"The CS 7200 fully meets my needs and workflow; it's very simple to use, with few buttons. Plus, it's small, compact and perfectly fits within the practice." Dr. Pascal Zyman



Plates that are as easy to position as film





High-dynamic range plates for reduced risk of under- or over-exposure



Same positioning and accessories as film

Available in most common film sizes



Exclusive Smart Plates



- Embedded RFID microchip identifies plates electronically
- Scan & Go workflow automates the scanning process
- Consistent quality and seamless performance with your Carestream Dental plate reader



Speeding up your workflow

FOUR EASY STEPS






Identify the plate electronically

2 Patient im the plate. Patient image is captured onto



Plate is scanned and erased 3 in one step. The plate is ready for reuse.








Scan & Go simplicity

SMART WORKFLOW FOR GREATER EFFICIENCY



- Truly designed for multiple patients and multiple users
- Multiple users can use the system at the same time without waiting in a queue
- Plates can be scanned in any order
- Eliminates risk of image mix ups and confusion

Fast and easy

- First image is ready to view in as little as eight seconds
- Scanner doesn't read if the plate is inserted upside down
- Plates are automatically erased after each scan
- Three scanning modes available to prioritize speed or resolution





Slim, space-saving design

ONLY 13 CM WIDTH SCANNER THAT FITS ANYWHERE








Fits into any dental practice

VERSATILE INSTALLATION



Compact design that can be used **chairside** with a direct USB connection



Unit can be **shared through a network** using Ethernet connection. No dedicated PC for the scanner needed



CS Imaging Version 8

Everything you need to see is now in just one place

- Streamline workflow—One-stop access to all your 2D images, 3D images and CAD/CAM data
- Facilitate collaboration—Easily share data with the treatment team
- Enhance diagnosis—Powerful analysis tools and cross-image comparison
- Increase case acceptance—Communicate treatment proposals more effectively
- Improve user experience—Modern and user-friendly design





CS Imaging Version 8

All your patient files in one place, instantly.

- Intraoral radiographs / photos
- Extraoral radiographs / photos
- CBCT volumes / screenshots
- Intraoral scanners files

All dental image formats and standards supported including JPEG, TIFF, DICOM, STL/PLY, TWAIN

As well as PDF, WORD, PPT, MP4






A more effective workflow

SEARCH, REVIEW AND SHARE IMAGES FROM ANYWHERE IN THE PRACTICE




Access to an instant overview of patient history and images



Select and enhance multiple images at the same time including FMX




Auto-arranges images so you don't have to manually move or resize images



Save and instantly retrieve all the images used during a treatment or the same visit



Quickly enhance clinical details using powerful processing tools




Drag and drop images to import and export (email, document or folder)



Confident diagnosis

NO COMPROMISE ON IMAGE QUALITY

- True resolution up to 19 lp/mm*
- Deliver consistently high-quality images
- Consistent quality with film-like rendering





Clinical results





Clinical results





The right accessories for the right results



Smart Plates Size 0, 1 or 2 package 4 plates per package



Hygienic sheaths Size 0, 1 or 2 package 200 pieces per package





Plate cleaning solution Plate s

Plate storage box and pad



Which system is right for you

QUICK COMPARISON






	CS 7200 Neo Edition	CS 7600
Use	Chairside / Shared use	Shared use / Multi-users
Theorical resolution (max.)	33 lp/mm	25 lp/mm
True resolution (max.)*	19 lp/mm	18 lp/mm
Smart Plates	~	~
Plate sizes	0-1-2	0-1-2-3-4
Examinations	Pediatric - Periapical - Bitewing	Pediatric - Periapical - Bitewing Long Bitewing - True occlusal
Time to first image	8 sec.	5 sec.
LCD screen with user instructions	-	√



Technical specifications

CS 7200	
Theoretical resolution (max)	33 lp/mm
True resolution (max)*	19 lp/mm
Imaging plate sizes (Height x width)	Size 0 – 22 mm x 35 mm Size 1 – 24 mm x 40 mm Size 2 – 31 mm x 41 mm
Plate erasing	Automatic
Scan time	As little as 8 seconds for one image
Connectivity	USB and Ethernet
Power supply	100–240 V (ac), 47 - 63 Hz, 0.8 - 0.4 A
System dimensions (W x H x D)	130 x 270 x 300 mm 5.1 x 10.6 x 11.8 in.
Weight	3.5 kg






Sales Information

- Product positioning
- Key selling points
- Competitive comparison
- Ordering information

Restricted to dealers and sales teams





Imaging plate systems portfolio

PRODUCT POSITIONING





Key Selling Points – CS 7200 Imaging Plate System

Advantages	Selling points	Proof points
High resolution images with a wider exposure range	 One of the highest true resolutions on the market (up to 19 lp/mm) High level of clinical details for confident diagnosis in any application, including endodontics Eliminates risks of under- or over-exposure 	
Easier and faster Scan & Go workflow	 Smart Plates are identified electronically prior to the exam and images are sent automatically to the right patient file Plates can be scanned in any order without any risk of plate mix-up Images are displayed in the correct position when scanning a Full Mouth Serie (FMS) 	
Fast scanning workflow	 First image is viewable in eight seconds and the plate automatically erased after reading Nearly instant results and diagnosis Limited waiting time between each scan 	
Easy to install and integrate	 The slim design (13 cm width) fits into tight spaces and can be installed chairside Easy-to-install system that can be installed by the user himself Flexible chairside or shared configuration that can adapt to any kind of practice–single operatory or multi-room 	
Reusable, film-like plates	 Thin, flexible plates are as easy to position as film and very comfortable for patients Scratch-resistant, durable plates that can be reused hundreds of times Plates available in same size as film 0,1 and 2 and covers most common indications 	
Powered with CS Imaging version 8	 One-stop access to all your images (2D, 3D, CAD/CAM) for a more effective workflow Easy image search, review and sharing from everywhere in the practice 	



Competitive comparison

Best	Good Poor or not ideal	and a second	de Dittore Vendeur	incon Trace		
	FEATURES	CS 7200 Neo Edition	Vistascan Nano Easy	Vistascan Mini Easy 2.0	Vistascan Mini view 2.0	Fussen
	Plate sizes	0,1,2	0,1,2	0,2	0,1,2,3,4	0,1,2,3
	Footprint	390cm2	360cm2	544cm2	544cm2	543cm2
	Connectivity	USB/Ethernet	Ethernet	Ethernet	Ethernet, WIFI	USB
	Smart workflow	Yes Patented technology	No	Yes	Yes	No
	Maximum resolution (theoretical)	33 lp/mm	16,7 lp/mm	40 lp/mm	40 lp/mm	12 lp/mm
	Strength /weakness				Al that orientate the exam automatically and monitors plate quality	



A REAL PROPERTY AND

442

Competitive comparison













FEATURES	CS 7200 Neo Edition	PSPix 2	Digora optime Dexis ScanExam	ScanX Duo	Apixia
Plate sizes	0,1,2	0,1,2,3 (+4 stich)	0,1,2,3	0,1,2,3,4	0,1,2,3
Footprint	390cm2	314cm2	699cm2	668cm2	252cm2
Connectivity	USB/Ethernet	Ethernet	Ethernet	USB	USB
Smart workflow	Yes Patented technology	No	No	No	No
Maximum resolution (theoretical)	33 lp/mm	20 lp/mm	17 lp/mm	40 lp/mm	17 lp/mm
Strength /weakness		Optical fiber to read the plate making sharper exams Smallest footprint of ht market			



CS 7200

DELIVERY CONTENT

Chairside configuration

- CS 7200 scanner
- Insertion panel size 0, 1, 2
- Plates size #2 (4 pieces)
- Hygienic sheaths size #2 (200 pieces)
- Screen Cleaner sample bottle (plate cleaning solution)
- USB cable
- Power supply and cord
- International AC adapter
- Quick user guide (poster)
- User and installation guide (USB key)

Network configuration

- CS 7200 scanner
- Scan & Go device (x2)
- Insertion panel size 0, 1, 2
- Plates size #2 (4 pieces)
- Hygienic sheaths size #2 (200 pieces)
- Screen Cleaner sample bottle (plate cleaning solution)
- Ethernet cable
- Power supply and cord
- International AC adapter
- Quick user guide (poster)
- User and installation guide (USB key)



Carestream's Dental imaging software should be ordered separately Additional accessories can be ordered separately



CAT Numbers

EQUIPMENT

	CAT #	Description	RRP
	5942727	CS 7200 IMAGING PLATE SYSTEM	€
NEW CAT#	5944608	CS 7200 + 2 x Scan & Go Devices	€

ACCESSORIES AND CONSUMABLES FOR CS 7200

\sim	CAT #	Description	RRP
	5942511	CSD Smart Plate #0 - 4 Pcs	€
	5942537	CSD Smart Plate #1 - 4 Pcs	€
	5942545	CSD Smart Plate #2 - 4 Pcs	€
	5942503	Hygienic Sheaths #0 - 200 Pcs	€
	5942529	Hygienic Sheaths #1 - 200 Pcs	€
	5942578	Hygienic Sheaths #2 - 200 Pcs	€
NEW CAT#	5944731	Scan & Go Device	
	1029594	CS 7600 / CS 7200 Transportation Box and Pad	€



Long-term peace of mind

- Two-year warranty
- Extended warranty option including parts replacement, offering you a long-term peace of mind

To be localized by regions



CS PROTECT

A warranty extension service to ensure peace of mind





Will the CS 7200 replace the CS 7600?

- No. Both products are positioned differently, with different features and different pricing
- The CS 7600 can read a broader range of plate sizes (including size 3 and 4) and features built-in memory and a color display.





Can we upgrade a CS 7200 reader to CS 7200 Neo Edition?

No. The former CS 7200 cannot be upgraded to a CS 7200 Neo Edition.

Can we adapt a CS 7200 Neo Edition with USB connection to a CS 7200 Neo Edition with Ethernet connection?

Yes. A CS 7200 Neo Edition with USB connection can be adapted with an Ethernet connection in order to share the reader through the network. This could be done by separately purchasing an ethernet cable (not provided by Carestream Dental) and Scan & Go devices to benefit from the full workflow.



Are the former CS 7200 plates and new plates interchangeable?

- The former CS 7200 plates are not compatible with the CS 7200 Neo Edition plate reader. Former plates cannot be read by a CS 7200 Neo Edition reader.
- The new CS 7200 / CS 7600 plates <u>are compatible</u> with the former CS 7200 plate reader but the Scan & Go capabilities won't be available when used with a former CS 7200 reader.



Former plates



New plates (for CS 7200 and CS 7600)



Are the new plates compatible with both the CS 7200 Neo Edition and CS 7600?

- Yes. The new plates are common to the CS 7200 Neo Edition and CS 7600 reader
- All the new Smart Plates feature an embedded Scan & Go identification chip that can be read by either CS 7200 Neo Edition and CS 7600 scanner







Are the hygienic sheaths for the CS 7200 Neo Edition similar to the ones for the CS 7600?

- No. The CS 7200 Neo Edition still uses sheaths that are specifically designed for the unit
- The sheaths for the two products open on different sides



CS 7200 sheaths



CS 7600 sheaths



How does the CS 7200 prevent cross-contamination?


- Several components prevent cross-contamination, including:
- Plates use single-use disposable hygienic sheaths
- Removable mechanical parts that are engaged with the imaging plate can be cleaned
- A new hygienic sheath is used for each use
- Cleaning and disinfecting plates between each patient is also recommended (cleaning solution)





Can we use the CS 7200 Neo Edition centrally?


- Yes. The CS 7200 Neo Edition is designed for both chairside use (paired with a dedicated workstation) or for shared used through the network.
- The main differences between two configurations are the connectivity (USB vs Ethernet) and the use of Scan & Go devices for the shared configuration.





What positioning accessories can be used with the Smart Plates?

- Plates can use the same positioning accessories as film
- Some manufacturers offer positioning devices that are better suited for imaging plates, such as plate holders designed to prevent damage to plates and to avoid punctures in hygienic sheaths
- "Alligator clips" should not be used due to risk of damaging plate and/or hygienic sheath






Can the CS 7200 unit be installed with software other than Carestream's Dental imaging software?

- Yes. TWAIN drivers are available for the CS 7200, which will allow users to review CS 7200 images within third-party, TWAIN-compliant software
- For more advanced software integration (SDK), please contact your Carestream Dental representative or the product marketing department












Scanner overview







Hygienic sheaths








Open a patient file in CS Imaging software.




Insert the plate into the sheath. Back side of plate should face opaque side of sheath. A B 1 Sec

Tag the plate with the Scan & Go device





Expose the plate



Open the sheath


Before insertion into the scanner, sheath should be opened by tearing the bottom side




6



Select plate size





Insert plate by:

- Placing the opaque side of the sheath up
- Aligning the left side of the sheath to the left side of the insertion slot





8



Pull the sheath to start the scanning process





Image appears in CS Imaging software



Thank You

